

Meicen A-Series AIO Baseplate

MC-AF001

- Provides multiple treatment options for head, neck and shoulders, thorax, pelvic and full-body support
- Different location choices to provide better immobilization for each patient
- Shoulder Retractor can be used to improve the accuracy and comfort during treatment

Specification: 1300(L) x 550(W)/mm



Also Available As MR Compatible

MC-AG001

MR



Meicen A-Series Head-Shoulder Baseplate

MC-AF002

- Provides immobilization solution for head, neck and shoulders, breast
- Sides emargination for simple installation and carry
- Contoured design for patient comfort
- Durable and lightweight

Specification: 1000(L) x 550(W)/mm



Also Available As MR Compatible

MC-AG002

MR




Meicen A-Series Baseplate

Pelvic Baseplate

MC-AF003

- Provides immobilization solution for thorax and abdomen
- Sides emargination for simple installation and carry
- Contoured design for patient comfort
- MR compatible MC-AG003

Specification: 850(L) x 520(W)/mm




Head-Shoulder Baseplate

MC-AF004

- Provides immobilization solution for head, neck and shoulders
- Compatible with all U-type Masks and L-type Masks
- Compatible with Color Headrest RE-HSP5
- MR compatible MC-AG004

Specification: 550(L) x 550(W)/mm



Compatible with



U-type Masks (Page 06-07)



Pelvic Masks



Headrests

Meicen W-Series Wing Board

MC-WBF001 T-Grip

MC-WBF002 U-Grip

- The Meicen Wing Board is durable and lightweight.
- It's easy to use with Vacuum Cushion greatly improves the comfort of the patient's treatment and the accuracy of positioning.
- The patient's arm is raised to maintain a grasping state, which can improve the stability and accuracy of the patient's immobilization.
- It can be used alone, or in conjunction with the Meicen C & Y-Series Baseplates, S-Series SBRT Solution and Supine Breast Board.



Wing Board T-Grip



Wing Board U-Grip

MC-WBF001-A01 T-Grip Carbon Fiber

MC-WBG001-A01 T-Grip MR



RC-HGF01 U-Grip Carbon Fiber

RM-HGG01 U-Grip MR



Also Available As MR Compatible

MC-WBG001

MR



Compatible with



T-Shaped Vacuum Cushion RC-E3010F