

Meicen B-Series


U-Frame Head Baseplate

MC-BF003

Carbon Fiber

- Provides the head support
- Compatible with all U-frame Masks
- Compatible with Meicen different material headrests
- MR compatible MC-BG003

Specification: 366(L) x 420(W)/mm




Tilting U-Frame Head Baseplate

MC-BF004

Carbon Fiber

- Provides the head support
- Compatible with all U-frame Masks
- Tilting angles: 5°/10°/20°/25°/30°/35°/40°/45°
- Compatible with Meicen different material headrests
- MR compatible MC-BG004

Specification: 377(L) x 306(W)/mm



Compatible with



U-type Masks (Page 34-35)



Indexing Bars