

Meicen R-Series **Breast Solutions**

MC-RBF001 Prone Breast Board

- Adjustable and durable massage-style face cushion.
- Multiple heterolateral wedge sizes to assist with patient tilt.
- Variable aperture openings to accommodate different breast sizes.
- Scale rulers for easy verification of longitudinal and lateral position of the patient.
- Reversible breast section allows for treatment of either breast.

Specification: 1570(L) x 450(W) x 150(H)/mm



Adjustable
Prone Headrest



MC-RBG001-A02



MC-RBG001-A03



MC-RBG001-A04



MC-RBG001-A05

Compatible with

M2422-VHT289W Trunk Mask 22% perforation 2.4mm

M3038-VHT550W Trunk Mask 38% perforation 3.0mm

Compatible with Trunk Mask together can be reduced the patient's involuntary displacement and improve the repeatability.

