

Meicen R-Series Breast Solutions

MC-RXE001 Supine Breast Cushions

The Meicen Supine Breast Cushions can be easily installed on Meicen C & Y-Series Baseplates by adapter cushion, or can also be indexed directly to couchtops. The Arm Support Cushions is ergonomically designed for patient comfort. The Hand Grips can be adjusted easily to improve the accuracy and comfort of repeatability.



Ergonomic design for patient comfort



Indexed slot for fixing on couchtops



Easy to install on Meicen Baseplates



The better positioning locations according to patient condition

Compatible with Meicen Baseplates



Meicen R-Series Breast Solutions



Wedge Cushion 8°
MC-RXE001-A01




Wedge Cushion 13°
MC-RXE001-A02



Arm Support Cushion (A-type)
MC-RXE001-A03



Arm Support Cushion (B-type)
MC-RXE001-A04



Arm Support Cushion (C-type)
MC-RXE001-A07



Adapter Cushion
MC-RXE001-A05



Color Headrest
RE-HSR5-B

Compatible with couchtops



Please contact our sales representative for more solutions.

Meicen R-Series Breast Solutions

MC-RBF001 Prone Breast Board



- Adjustable and durable massage-style face cushion.
- Multiple heterolateral wedge sizes to assist with patient tilt.
- Variable aperture openings to accommodate different breast sizes.
- Scale rulers for easy verification of longitudinal and lateral position of the patient.
- Reversible breast section allows for treatment of either breast.

Specification: 1570(L) x 450(W) x 150(H)/mm



MC-RBG001-A02



MC-RBG001-A03



MC-RBG001-A04





MC-RBG001-A05

Compatible with

M2422-VHT289W Trunk Mask 22% perforation 2.4mm

M3038-VHT550W Trunk Mask 38% perforation 3.0mm

Compatible with Trunk Mask together can be reduced the patient's involuntary displacement and improve the repeatability.



Meicen R-Series Breast Solutions

MC-RXF002 Supine Breast Board



- 5 inclinations: 5°, 10°, 15°, 20° and 25°.
- Simple, customizable setups for patient comfort and clearance.
- Contoured lumbar area and enhanced bottomstop with easy locking.
- Compatible with S-type Head Masks for stable, reproducible supraclavicular treatments.

Specification: 1325(L) x 445(W)/mm

Also Available As MR Compatible

MC-RXG002

MR



Wrist Support
MC-RXF001-A01
MC-RXG001-A01



Arm Support
MC-RXF001-A02
MC-RXG001-A02

Compatible with

M2422-VSH190W Breast Mask 22% perforation 2.4mm

M3038-VSH480W Breast Mask 38% perforation 3.0mm

Compatible with Breast Mask can be reduced the patient's involuntary displacement and improve the repeatability and comfort of the positioning.



Meicen R-Series Breast Solutions

MC-RBE001 Prone Breast Cushion

- The Breast-bridge can be converted left and right according to the treatment plan.
- The U-shaped cushion can be selected S, M, L size.
- The rail-type prone headrest with head to feet direction and the four-levels elevation angle can be adjusted to meet individual needs.
- Compatible with Meicen C & Y-Series AIO Baseplates.



Compatible with

M2422-VST190W Trunk Mask 22% perforation 2.4mm

M3038-VST480W Trunk Mask 38% perforation 3.0mm

Compatible with Trunk Mask together can be reduced the patient's involuntary displacement and improve the repeatability.



Meicen R-Series Prone Belly Solutions

MC-RPE002 Prone Belly Cushion

- The rail-type prone headrest with head to feet direction and the four-levels elevation angle can be adjusted to meet individual needs.
- Dual-configured for female and male treatments.
- Compatible with Meicen C & Y-Series AIO Baseplates.



Prone Headrest



Compatible with

M2422-VST190W Trunk Mask 22% perforation 2.4mm

M3038-VST480W Trunk Mask 38% perforation 3.0mm

Compatible with Trunk Mask together can be reduced the patient's involuntary displacement and improve the repeatability.

