

Meicen Vacuum Cushions

The Meicen Vacuum Cushions offer a rigid and secure support around the patients. The coated nylon shell protects them from tear and stain. The Vacuum Cushions are completely air-tight, leak proof, easy clean and fill with optional density particles made from polyamide. Built-in indexing bars are available, different size and shape of Vacuum Cushions are available too.



- Ultra durable, high-quality polyamide
- Radio-translucent
- Easy to clean
- Cost-effective
- Adapters available for any manufacturer's pump
- Dark blue color distinguishable in optical imagers

Meicen Vacuum Cushions




Non-indexed



Indexed



CT Valve



MR Valve



Indexing Bar



Adapter

Rectangular Vacuum Cushions

CT	MR	Dimensions	Fill	Indexed	
RC-E5003 RC-E5005S	RM-E5003 RM-E5005S	100x70cm	30L	W/Indexed	for chest-pelvic support
RC-E4013 RC-E4014F	RM-E4013 RM-E4014F	70x70cm	20L	W/Indexed	for breast support
RC-E5024 RC-E5025S	RM-E5024 RM-E5025S	120x80cm	32L	W/Indexed	for chest-pelvic support
RC-E5010 RC-E5012S	RM-E5010 RM-E5012S	100x75cm	30L	W/Indexed	for chest-pelvic support
RC-E5023S	RM-E5023S	120x75cm	34L	W/Indexed	for chest-pelvic support
RC-E5032 RC-E5033S	RM-E5032 RM-E5033S	150x75cm	45L	W/Indexed	for chest-pelvic support
RC-E8005	RM-E8005	80x200cm	80L		for full-body support
RC-E8020S	RM-E8020S	100x150cm	75L	W/Indexed	for full-body support

Remark: Custom width shall not exceed 140cm, Custom length shall not exceed 250cm.

Meicen Vacuum Cushions



T-Shape Vacuum Cushion
with Open Headrest
RC-E3010F



T-Shape Vacuum Cushion
RC-E3005F



Wing Board Cushion
RC-E4003



S-Shape Vacuum Cushion
RC-E1004

Non-Rectangular Vacuum Cushions

CT	MR	Dimensions	Fill	Shape	
RC-E1005	RM-E1005	75x75cm	15L	Arms up	for breast support
RC-E3005F RC-E3010F	RM-E3005F RM-E3010F	125x110cm	31L 32L	T-Shaped	for thorax-abdomen support
RC-E3019	RM-E3019	105x75cm	28L	Arms up	for thorax-abdomen support
RC-E4003 RC-E4004	RM-E4003 RM-E4004	110x55cm	15L 19L	Wing Board breast	for breast support
RC-E1004	RM-E1004	61x65cm	8L	S-Shaped	for head-shoulder support
RC-E1011F	RM-E1011F	68x61cm	10L	S-Shaped	for head-shoulder support
RC-E1008	RM-E1008	50x70cm	12L	S-Shaped	for head-shoulder support