

## Meicen Y-Series AIO Baseplate

The Meicen Y-Series Baseplate have a Pin-lock design for easy operation, provide the highest degree of stability for brain, head and neck, thorax, pelvic and full-body support. Dual-configured for adult and pediatric treatments. It is compatible with all S-type Masks and Meicen Pin-lock type Masks.



Y-Series AIO Baseplate Solutions

### Optional Accessories



Carbon Fiber Headrest  
RE-HSF6



Overhead Arm Support  
RC-ASF01



Shoulder Retractor  
MC-YG001-A03




Belly Bridge  
MC-YF001-A01

## Meicen Y-Series AIO Baseplate

### MC-YF001

- Compatible with standard and narrow shoulder S-type Masks
- Combined with Belly Bridge to form a SBRT solution
- Different location choices to provide better immobilization for each patient
- Configured for multiple treatments


Specification: 1300(L) x 606(W)/mm



Pins for fixing headrests



Easily adjustable Hand Grips



Adjustable to patient size



Optional locations for Leg Separator  
MC-CG001-A04

Also Available As MR Compatible

### MC-YG001



### Compatible with



Indexing Bars

## Meicen Y-Series AIO Baseplate

### Optional Accessories

#### Belly Bridge MC-YF001-A01

- Ergonomic Compression Paddle for better comfort.
- The paddle device can be adjusted in the direction of left and right, up and down, and has a scale to record its fixed position.
- Abundant scale marks improve the accuracy of positioning repeatability.
- MR compatible MC-YG001-A01



#### Knee Bridge MC-YF001-A02

- The knee paddle device can be adjusted left and right, up and down, and has scales to record its fixed position.
- Abundant scale marks improve the accuracy of positioning repeatability.
- Knee paddle device is optional in both legs or single leg immobilization.
- MR compatible MC-YG001-A02



#### Extension Board RC-EBF01

- Improve the comfort of the patient when positioning.
- Can be fixed directly onto the couchtop.
- MR compatible RM-EBG01



#### Overhead Arm Support RC-ASF01

- It can be used together with the Meicen C & Y-Series Baseplates to improve the reproducible precision and comfort in breast thorax and lung tumor treatment.
- The height of Wrist Supports and angle of Arm Supports can be adjusted to each patient.
- MR compatible RM-ASG01



Wrist Support  
RC-ASF01-A02



Arm Support  
RC-ASF01-A01

## Meicen Y-Series AIO Baseplate



### SBRT Solution Package

Components	CF Package	MR Package	Qty
Belly Bridge	MC-YF001-A01	MC-YG001-A01	1
Knee Bridge	MC-YF001-A02	MC-YG001-A02	1
Extension Board	RC-EBF01	RM-EBG01	1
Overhead Arm Support	RC-ASF01	RM-ASG01	1
Knee Cushion	RE-FCA01	RE-FCA01	1
Elevation Cushion	RE-FCA02	RE-FCA02	1
Feet Cushion	RE-FCA03	RE-FCA03	1
Headrest	RE-HSF6	RE-HSB5	1
Indexing Bar (for E-style indexed couchtops)	RC-LBE02	RM-LBE01	3
Indexing Bar (for V-style indexed couchtops)	RC-LBV02	RM-LBV01	3
Shoulder Retractor	MC-YG001-A03	MC-YG001-A03	1
SBRT Indexing Bar	MC-STF001-A07	MC-STG001-A07	2

### Also Available As MR Compatible

