

## Meicen Violet Thermoplastic Masks


Meicen Violet Thermoplastic Masks are composed of a new, premium thermoplastic material designed specifically to offer the ultimate in patient comfort with exceptional rigidity and extremely low shrinkage. Available in multiple perforation styles and compatible with most manufacturers' baseplates, our violet thermoplastic's primary focus is the patient experience.

- Low shrinkage
- Optimal molding time
- Ergonomically designed for patient comfort
- Non-Stick
- Violet color to alleviate patient anxiety




## Meicen Violet Thermoplastic Masks


Meicen Reinforced Series Masks feature a combination of enhanced perforation and solid thermoplastic to provide rigidity and secure immobilization for IMRT treatments. Openface Masks can be compatible with SGRT (Surface Guided Radiation Therapy) camera systems.




S-type Reinforced



P-type Reinforced



L-type Reinforced



Stereotactic X-Knife

Meicen Violet Thermoplastic Masks are easy to prepare and can be heated by water bath at 65-70°C (149-158°F) or oven at 70-75°C (158-167°F).

Moulding Time: 1-2 minutes



Cooling Time: 8-15 minutes

Thickness 2.0mm, 2.4mm, 3.0mm or 3.2mm are available.

Please contact our sales representative for more information.

## U-frame Violet Thermoplastic Masks

- Slim Frame



### U-Shaped Head Mask

Standard


M2422-VUA192W	22% perforation, 2.4mm
M3038-VUA481W	38% perforation, 3.0mm
M2415-VUA150W	15% perforation, 2.4mm



### U-Shaped Head Mask

Reinforced


M24RA-VUA580W	2.4mm
M30RA-VUA610W	3.0mm



### U-Shaped Head Mask

Reinforced, Openface

M24RB-VUA640W	2.4mm
M30RB-VUA670W	3.0mm



### U-Shaped Head Mask

Reinforced, Open eyes & mouth

M24RB-VUA641W	2.4mm
M30RB-VUA671W	3.0mm




### U-Shaped Head Mask

Standard, Extended

		Extended 5cm	Extended 9cm	
M2422-VUA197W	M2422-VUA200W			22% perforation, 2.4mm
M3038-VUA482W	M3038-VUA483W			38% perforation, 3.0mm

## U-frame Violet Thermoplastic Masks

- Standard Frame



### U-Shaped Head Mask

Standard


M2422-VUA190W	22% perforation, 2.4mm
M3038-VUA480W	38% perforation, 3.0mm



### U-Shaped Head Mask

Reinforced


M24RA-VUA581W	2.4mm
M30RA-VUA611W	3.0mm



### U-Shaped Head Mask

Reinforced, Openface

M24RB-VUA642W	2.4mm
M30RB-VUA672W	3.0mm



### U-Shaped Head Mask

Reinforced, Open eyes & mouth

M24RB-VUA643W	2.4mm
M30RB-VUA673W	3.0mm



### U-Shaped Head Mask

Standard, Extended


		Extended 5cm	Extended 9cm	
M2422-VUA194W	M2422-VUA202W			22% perforation, 2.4mm
M3038-VUA484W	M3038-VUA485W			38% perforation, 3.0mm

Meicen U-frame Masks fit all standard U-frame baseplates from Meicen and other manufacturers.

2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.




## ► S-type Violet Thermoplastic Masks



### S-Shaped Head Mask

Standard


M2422-VSA190W	22% perforation, 2.4mm
M3038-VSA480W	38% perforation, 3.0mm
M2415-VSA150W	15% perforation, 2.4mm



### S-Shaped Head Mask

Reinforced

M24RA-VSA580W	2.4mm
M30RA-VSA610W	3.0mm



### S-Shaped Head Mask

Standard, Extended

Extended 5cm		Extended 9cm		
M2422-VSA192W	M2422-VSA193W	22% perforation, 2.4mm		
M3038-VSA482W	M3038-VSA483W	38% perforation, 3.0mm		




### S-Shaped Head Mask

Reinforced, Extended

Extended 5cm		Extended 9cm		
M24RA-VSA581W	M24RA-VSA582W	2.4mm		
M30RA-VSA611W	M30RA-VSA612W	3.0mm		

Meicen S-type Masks fit all standard S-type baseplates from Meicen and other manufacturers.  
2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.


## ► S-type Violet Thermoplastic Masks



### S-Shaped Head-Neck-Shoulder Mask

Standard


M2422-VSB190W	22% perforation, 2.4mm
M3038-VSB480W	38% perforation, 3.0mm
M2415-VSB150W	15% perforation, 2.4mm



### S-Shaped Head-Neck-Shoulder Mask

Reinforced


M24RA-VSB580W	2.4mm
M30RA-VSB610W	3.0mm



### S-Shaped Head-Neck-Shoulder Mask (narrow shoulders)

Standard

M2422-VSB191W	22% perforation, 2.4mm
M3038-VSB481W	38% perforation, 3.0mm



### S-Shaped Head-Neck-Shoulder Mask (narrow shoulders)

Reinforced

M24RA-VSB581W	2.4mm
M30RA-VSB611W	3.0mm


\* Compatible with Meicen C & Y-Series baseplates.

## ► S-type Violet Thermoplastic Masks (Openface)




Meicen Openface Series Masks aims to alleviate symptoms of claustrophobia during radiotherapy. All types of masks are available as openfaced.

Adjustable pulling of the open part when apply to different patients assures comfort and reinforced effect.




**S-Shaped Head Mask**  
Reinforced, Open eyes & mouth

M24RB-VSA643W	2.4mm
M30RB-VSA672W	3.0mm



**S-Shaped Head Mask**  
Reinforced, Openface


M24RB-VSA640W	2.4mm
M30RB-VSA670W	3.0mm



**S-Shaped Head Mask**  
Reinforced, Openface


M24RB-VSA641W	2.4mm
M30RB-VSA671W	3.0mm

## ► S-type Violet Thermoplastic Masks (Openface)




**S-Shaped Head-Neck-Shoulder Mask**  
Reinforced, Openface

M24RB-VSB642W	2.4mm
M30RB-VSB670W	3.0mm




**S-Shaped Head-Neck-Shoulder Mask**  
Reinforced, Open eyes & mouth

M24RB-VSB643W	2.4mm
M30RB-VSB672W	3.0mm



**S-Shaped Head-Neck-Shoulder Mask**  
Reinforced, Openface

M24RB-VSB641W	2.4mm
M30RB-VSB671W	3.0mm




**S-Shaped Head-Neck-Shoulder Mask**  
Reinforced, Open face & neck

M24RB-VSB644W	2.4mm
M30RB-VSB673W	3.0mm

Meicen S-type Masks fit all standard S-type baseplates from Meicen and other manufacturers.  
2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.



## ► L-type Violet Thermoplastic Masks



### Triangular Head Mask

Standard, with nose hole

M2422-VLA191W	22% perforation, 2.4mm
M3038-VLA480W	38% perforation, 3.0mm
M2438-VLA320W	38% perforation, 2.4mm, for pediatric



### Triangular Head Mask

Standard, Extended 3cm, with nose hole


M2422-VLA192W	22% perforation, 2.4mm
M3038-VLA481W	38% perforation, 3.0mm
M2415-VLA151W	15% perforation, 2.4mm



### Triangular Head Mask

Reinforced, with nose hole


M24RA-VLA580W	2.4mm
M30RA-VLA610W	3.0mm



### Neck-Shoulder Mask

Standard, with nose hole

M2422-VLE190W	22% perforation, 2.4mm
M3038-VLE480W	38% perforation, 3.0mm
M2422-VLE191W	22% perforation, 2.4mm, for pediatric




### Neck-Thorax Mask

Standard

M2422-VLF192W	22% perforation, 2.4mm
M3038-VLF480W	38% perforation, 3.0mm


## ► L-type Violet Thermoplastic Masks



### Head-Neck-Shoulder Mask

Standard, with nose hole


M2422-VLB191W	22% perforation, 2.4mm
M3038-VLB481W	38% perforation, 3.0mm
M2422-VLB192W	22% perforation, 2.4mm, for pediatric



### Head-Neck-Shoulder Mask

Reinforced, with nose hole


M24RA-VLB580W	2.4mm
M30RA-VLB610W	3.0mm



### Chest-Pelvic Mask

Standard


M3038-VLK483W	38% perforation, 3.0mm, 228 clamp
M3038-VLK482W	38% perforation, 3.0mm, 243 clamp
M3038-VLK481W	38% perforation, 3.0mm, 247 clamp



### Pelvic Mask

Standard

M3038-VLL483W	38% perforation, 3.0mm, 228 clamp
M3038-VLL482W	38% perforation, 3.0mm, 243 clamp
M3038-VLL481W	38% perforation, 3.0mm, 247 clamp




### Thorax-Abdomen Mask

Standard

M2422-VLG190W	22% perforation, 2.4mm
M3038-VLG480W	38% perforation, 3.0mm


Meicen L-type Masks fit all standard O-type baseplates from Meicen and other manufacturers.  
2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.

## ► L-type Violet Thermoplastic Masks (Openface)




**Triangular Head Mask**  
Reinforced, Openface

M24RB-VLA641W	2.4mm
M30RB-VLA671W	3.0mm




**Triangular Head Mask**  
Reinforced, Open eyes & mouth

M24RB-VLA643W	2.4mm
M30RB-VLA673W	3.0mm




**Head-Neck-Shoulder Mask**  
Reinforced, Open eyes & mouth

M24RB-VLB642W	2.4mm
M30RB-VLB672W	3.0mm



**Head-Neck-Shoulder Mask**  
Openface

M24RB-VLB643W	2.4mm
M30RB-VLB673W	3.0mm



**Head-Neck-Shoulder Mask**  
Reinforced, Extended, Openface

M24RB-VLC641W	2.4mm
M30RB-VLC671W	3.0mm

Meicen L-type Masks fit all standard O-type baseplates from Meicen and other manufacturers.  
2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.

## ► P-type Violet Thermoplastic Masks




**Triangular Head Mask**  
Standard, with nose hole

M2422-VPA190W	22% perforation, 2.4mm
M3038-VPA480W	38% perforation, 3.0mm
M2415-VPA150W	15% perforation, 2.4mm




**Triangular Head Mask**  
Reinforced, with nose hole

M24RA-VPA580W	2.4mm
M30RA-VPA610W	3.0mm




**Triangular Head Mask**  
Standard, Extended 5cm, with nose hole

M2422-VPA191W	22% perforation, 2.4mm
M3038-VPA481W	38% perforation, 3.0mm
M2415-VPA151W	15% perforation, 2.4mm



**Neck-Shoulder Mask**  
Standard, with nose hole


M2422-VPE190W	22% perforation, 2.4mm
M3038-VPE480W	38% perforation, 3.0mm



**Head-Neck-Shoulder Mask**  
Standard, with nose hole

M2422-VPB190W	22% perforation, 2.4mm
M3038-VPB480W	38% perforation, 3.0mm
M2415-VPB150W	15% perforation, 2.4mm

## ► P-type Violet Thermoplastic Masks



### Head-Neck-Shoulder Mask

Reinforced, with nose hole

M24RA-VPB580W	2.4mm
M30RA-VPB610W	3.0mm



### Chest-Pelvic Mask

Standard

M2422-VPK190W	22% perforation, 2.4mm
M3038-VPK481W	38% perforation, 3.0mm



### Triangular Head Mask (for pediatric)

Standard, with nose hole


M2422-VPA192W	22% perforation, 2.4mm
M3038-VPA482W	38% perforation, 3.0mm



### Head-Neck-Shoulder Mask (for pediatric)

Standard, with nose hole

M2422-VPB191W	22% perforation, 2.4mm
M3038-VPB481W	38% perforation, 3.0mm



### Head-Neck-Shoulder Mask (for pediatric)

Reinforced, with nose hole

M24RA-VPB581W	2.4mm
M30RA-VPB611W	3.0mm

## ► P-type Violet Thermoplastic Masks (Openface)



### Triangular Head Mask

Standard, Openface

M24RB-VPA640W	2.4mm
M30RB-VPA670W	3.0mm



### Triangular Head Mask

Reinforced, Openface


M24RB-VPA641W	2.4mm
M30RB-VPA671W	3.0mm



### Triangular Head Mask

Reinforced, Open eyes & mouth


M24RB-VPA642W	2.4mm
M30RB-VPA672W	3.0mm



### Head-Neck-Shoulder Mask

Reinforced, Openface

M24RB-VPB640W	2.4mm
M30RB-VPB670W	3.0mm



### Head-Neck-Shoulder Mask

Reinforced, Open eyes & mouth

M24RB-VPB641W	2.4mm
M30RB-VPB671W	3.0mm


Meicen P-type Masks fit all PosiFix-type baseplates from Meicen and Civco.

2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.



## Pin-lock type Violet Thermoplastic Masks


- For the Meicen C & Y-Series baseplates.



### Neck-Thorax Mask

Reinforced


M24RA-VSF580W	2.4mm
M30RA-VSF610W	3.0mm



### Trunk Mask (Chest-Pelvic)

Standard

M2422-VST190W	22% perforation, 2.4mm
M3038-VST480W	38% perforation, 3.0mm




### Trunk Mask (Pelvic)

Standard, with Groin Lock

M2422-VST191W	22% perforation, 2.4mm
M3038-VST481W	38% perforation, 3.0mm

For Meicen Supine Breast Board



### Breast Mask


Standard

M2422-VSH190W	22% perforation, 2.4mm
M3038-VSH480W	38% perforation, 3.0mm

2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.

## H-type Violet Thermoplastic Masks


- For the Civco HipFix System.



### Trunk Mask (Chest-Pelvic)

Standard


M2422-VHT190W	22% perforation, 2.4mm
M3038-VHT480W	38% perforation, 3.0mm



### Trunk Mask (Pelvic)

Standard, with Groin Lock

M2422-VHT191W	22% perforation, 2.4mm
M3038-VHT481W	38% perforation, 3.0mm



### Trunk Mask (Chest-Pelvic)

Reinforced

M24RA-VHT580W	2.4mm
M30RA-VHT610W	3.0mm



### Trunk Mask (Chest-Pelvic)


Plain

M2400-VHT120W	2.4mm
M3000-VHT420W	3.0mm

2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.

## M-type Violet Thermoplastic Masks

- For the MacroMedics OmniBoard 2.0



### Triangular Head Mask

Standard, with nose hole


M2422-VMA190W	22% perforation, 2.4mm
M3038-VMA480W	38% perforation, 3.0mm



### Triangular Head Mask

Reinforced, with nose hole


M24RA-VMA580W	2.4mm
M30RA-VMA610W	3.0mm



### Head-Neck-Shoulder Mask

Standard, with nose hole


M2422-VMB190W	22% perforation, 2.4mm
M3038-VMB480W	38% perforation, 3.0mm



### Head-Neck-Shoulder Mask

Reinforced, with nose hole

M24RA-VMB580W	2.4mm
M30RA-VMB610W	3.0mm



### Neck-Shoulder Mask

Standard


M2422-VME190W	22% perforation, 2.4mm
M3038-VME480W	38% perforation, 3.0mm

2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.

## X-Knife Violet Thermoplastic Masks


Meicen X-Knife Violet Thermoplastic Masks are composed of a new, premium thermoplastic material designed specifically to offer the ultimate in patient comfort with exceptional rigidity and extremely low shrinkage. Additionally, the X-Knife Violet Thermoplastic Masks have a non-stick coating applied to the patient-facing mask surfaces, helping ensure that the masks do not stick to the patient's hair, beard, etc.

- Adhesive Layers
- Low shrinkage
- Easy to manipulate and mould
- Cost-effective




### X-Knife Mask

38% perforation, 3.0mm  
M3038-VXA481W



### X-Knife Head-Shoulder Mask

38% perforation, 3.0mm  
M3038-VXB480W




### X-Knife Mask

38% perforation, 3.0mm  
M3038-VXA490W

## E-type Violet Thermoplastic Masks


Compatible with Elekta Fraxion System



### Triangular Head Mask

Reinforced

M24RA-VEA600W	2.4mm
M30RA-VEA620W	3.0mm



### Triangular Head Mask

Reinforced

M24RA-VEA602W	2.4mm
M30RA-VEA621W	3.0mm

Compatible with Elekta HeadSTEP & HipSTEP Baseplates



### Head Mask

Standard

M2422-VEA190W	22% perforation, 2.4mm
M3038-VEA480W	38% perforation, 3.0mm
M2415-VEA150W	15% perforation, 2.4mm




### Head Mask

Reinforced

M24RA-VEA580W	2.4mm
M30RA-VEA610W	3.0mm

2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.

## E-type Violet Thermoplastic Masks



### Head Mask

Standard, Extended 5cm


M2422-VEA192W	22% perforation, 2.4mm
M3038-VEA482W	38% perforation, 3.0mm



### Head Mask

Standard, Extended 21cm


M2422-VEA193W	22% perforation, 2.4mm
M3038-VEA483W	38% perforation, 3.0mm



### Head-Neck-Shoulder Mask

Standard

M2422-VEB190W	22% perforation, 2.4mm
M3038-VEB480W	38% perforation, 3.0mm
M2415-VEB150W	15% perforation, 2.4mm



### Head-Neck-Shoulder Mask

Reinforced

M24RA-VEB580W	2.4mm
M30RA-VEB610W	3.0mm



### Chest-Pelvic Mask

Standard

M3038-VEK480W	38% perforation, 3.0mm
---------------	------------------------




## ➤ V-type Violet Thermoplastic Masks



### V-Shaped Head Mask

Standard


M2422-VVA190W	22% perforation, 2.4mm
M3038-VVA480W	38% perforation, 3.0mm



### V-Shaped Head Mask

Reinforced


M24RA-VVA580W	2.4mm
M30RA-VVA610W	3.0mm



### V-Shaped Head Mask

Reinforced, Open eyes & mouth


M24RB-VVA640W	2.4mm
M30RB-VVA670W	3.0mm



### V-Shaped Head Mask

Reinforced, Openface

M24RB-VVA641W	2.4mm
M30RB-VVA671W	3.0mm




### V-Shaped Head Mask

Standard, Extended

		Extended 5cm	Extended 9cm	
M2422-VVA191W	M2422-VVA192W			22% perforation, 2.4mm
M3038-VVA481W	M3038-VVA482W			38% perforation, 3.0mm


## ➤ V-type Violet Thermoplastic Masks



### V-Shaped Head-Neck-Shoulder Mask

Standard


M2422-VVB190W	22% perforation, 2.4mm
M3038-VVB480W	38% perforation, 3.0mm



### V-Shaped Head-Neck-Shoulder Mask

Reinforced


M24RA-VVB580W	2.4mm
M30RA-VVB610W	3.0mm



### V-Shaped Head-Neck-Shoulder Mask

Reinforced, Openface

M24RB-VVB640W	2.4mm
M30RB-VVB670W	3.0mm



### V-Shaped Head-Neck-Shoulder Mask

Reinforced, Open eyes & mouth

M24RB-VVB641W	2.4mm
M30RB-VVB671W	3.0mm

For the Klarity Optek board and Bionix VersaBoard.

2.0mm & 3.2mm thickness are available also, please specify thickness when requesting a quote.