

Meicen R-Series Breast Solutions

MC-RXF002 Supine Breast Board

- 5 inclinations: 5°, 10°, 15°, 20° and 25°.
- Simple, customizable setups for patient comfort and clearance.
- Contoured lumbar area and enhanced bottomstop with easy locking.
- Compatible with Meicen Breast Masks.
- Compatible with S-type Head Masks for stable, reproducible supraclavicular treatments.

Specification: 1325(L) x 445(W)/mm



Also Available As MR Compatible

MC-RXG002

MR



Wrist Support
MC-RXF001-A01
MC-RXG001-A01



Arm Support
MC-RXF001-A02
MC-RXG001-A02

Compatible with

M2422-VSH190W Breast Mask 22% perforation 2.4mm

M3038-VSH480W Breast Mask 38% perforation 3.0mm

Compatible with Breast Mask can be reduced the patient's involuntary displacement and improve the repeatability and comfort of the positioning.

