

Meicen Vacuum Cushions

The Meicen Vacuum Cushions offer a rigid and secure support around the patients. The coated nylon shell protects them from tear and stain. The Vacuum Cushions are completely air-tight, leak proof, easy clean and fill with optional density particles made from polyamide. Built-in indexing bars are available, different size and shape of Vacuum Cushions are available too.



Meicen Vacuum Cushions

- Ultra durable, high-quality polyamide.
- Dark / Light blue color distinguishable in optical imagers.
- Radio-translucent.
- Easy to clean.
- Cost-effective.
- Adapters available for any manufacturer's pump.



Valve



Indexing Bar



CPC Adapter

CT	MR	Dimensions	Fill	Indexed	
RC-E5003 RC-E5005S	RM-E5003 RM-E5005S	100x70cm	30L	W/Indexed	For Chest-Pelvic Support
RC-E4013 RC-E4014F	RM-E4013 RM-E4014F	70x70cm	20L	W/Indexed	For Breast Support
RC-E5024 RC-E5025S	RM-E5024 RM-E5025S	120x80cm	32L	W/Indexed	For Chest-Pelvic Support
RC-E5010 RC-E5012S	RM-E5010 RM-E5012S	100x75cm	30L	W/Indexed	For Chest-Pelvic Support
RC-E5023S	RM-E5023S	120x75cm	34L	W/Indexed	For Chest-Pelvic Support
RC-E5032 RC-E5033S	RM-E5032 RM-E5033S	150x75cm	45L	W/Indexed	For Chest-Pelvic Support

Remark: Custom width shall not exceed 140cm

Custom length shall not exceed 250cm

please specify color (Dark blue or Light blue color) when requesting a quote.

Meicen Vacuum Cushions



T-Shape Vacuum Cushion
with Open Headrest, 125 x 110cm
RC-E3010F



T-Shape Vacuum Cushion
125 x 110cm
RC-E3005F



Wing Board Cushion
110 x 55cm
RC-E4003



S-Shape Vacuum Cushion
61 x 65cm
RC-E1004

Non-Rectangular Vacuum Cushions

CT	MR	Dimensions	Fill	Shape	
RC-E1005	RM-E1005	75x75cm	15L	Arms up	For Breast Support
RC-E3005F RC-E3010F	RM-E3005F RM-E3010F	125x110cm	31L 32L	T-Shaped	For Thorax-Abdomen Support
RC-E3019	RM-E3019	105x75cm	28L	Arms up	For Thorax-Abdomen Support
RC-E4003 RC-E4004	RM-E4003 RM-E4004	110x55cm	15L 19L	Wing board breast	For Breast Support
RC-E1004	RM-E1004	61x65cm	8L	S-Shaped	For Head-Shoulder Support

Diagnostic Couchtop Overlays

- The Couchtop Overlays provide a clean, flat surface for accurate and easily reproducible patient positioning.
- The internal density distribution of the product is uniform, and there is no unacceptable artifact in the scanning area.
- Can be quickly installed and disassembled, simple and convenient.
- It is stable and reliable after installation and meets the leveling requirements.
- The Couchtop Overlays are designed to work with Indexing bars as well as other positioning systems.


* The Couchtop Overlays can be customized according to the user's CT or MRI equipment.

MC-CTF019 CT Overlay



Compatible to Somatom CT Definition AS+

MC-MRG002 MRI Overlay



Compatible to Philips Ingenia Overlay